# Occlutech

INSTRUCTIONS FOR USE

**Figulla® Flex II ASD Occluder**

<b>English</b> .....	<b>3</b>
<b>Deutsch</b> .....	<b>17</b>
<b>Español</b> .....	<b>31</b>
<b>French</b> .....	<b>45</b>
<b>Italiano</b> .....	<b>59</b>
<b>Polski</b> .....	<b>73</b>
<b>Português</b> .....	<b>88</b>
<b>Türkçe</b> .....	<b>102</b>

# TABLE OF CONTENTS


1. Product description .....	3
2. Indication for Use and Area of Application .....	4
3. Contraindications.....	4
4. Warnings .....	5
5. Preventative measures .....	6
5.1 Handling .....	6
5.2 Device Size.....	7
6. Adverse events.....	7
7. Magnetic Resonance Imaging (MRI) safety information .....	7
8. Implantation Procedure.....	8
8.1 Measurement of Defect and Device Sizing .....	8
8.2 Procedural Method of Implantation .....	8
8.3 Notes for Troubleshooting .....	12
8.4 Report of Serious Adverse Events.....	13
8.5 Recommended Follow-up Examinations.....	13
8.6 Recommended Medication after Implantation .....	13
9. Device Sizes and Recommended delivery systems .....	14
10. Storage conditions .....	15
11. Disposal .....	15
12. Labelling.....	16

---

## 1. PRODUCT DESCRIPTION

The Figulla® Flex II Atrial Septal Defect (ASD) Occluder is a percutaneous transcatheter closure device consisting of two Nitinol wire mesh discs connected by a flexible waist. The design properties allow the device to be pulled into a delivery sheath and reshape once it is released. When the device is correct positioned and firmly attached on both sides of the septum, it conforms to the septal wall and closes a present ASD. Depending on the Figulla® Flex II ASD Occluder size, two or three very thin polyester patches, sewn into the wire mesh, are intended to stop the blood flow through the ASD as well as optimize tissue growth.

The Figulla® Flex II ASD Occluder is packed together with its compatible pusher as shown in Figure 1.

**Figure 1:** (A) – Product box as delivered to the physician; (B) – Figulla® Flex II ASD Occluder; (C) – Patient Card; (D) – pusher (for example: Occlutech® Pistol Pusher); (E) – Instructions for Use (one for the Figulla® Flex II ASD Occluder and one for the Occlutech® Pistol Pusher).

## 2. INDICATION FOR USE AND AREA OF APPLICATION

The Figulla® Flex II ASD Occluder is a medical device intended for transcatheter closure of ostium secundum-type ASD. Patients indicated for ASD closure have:

- Echocardiographic evidence of ostium secundum-type ASD
- Clinical evidence of right ventricular (RV) volume overload (hemodynamically significant left-to-right shunt with  $Q_p / Q_s \geq 1.5$  or RV enlargement)

## 3. CONTRAINDICATIONS

The Figulla® Flex II ASD Occluder is contraindicated for the following:

- Acute infection
- Adequate oral anti-coagulation therapy / platelet inhibition is not possible post-intervention
- Arrhythmia
- ASD larger than 40 mm
- ASD primum defect
- Atrial tumor
- Atrial-thrombus
- Eisenmenger-syndrome
- Known coagulation disorder
- Intolerance to contrast agents
- Allergies to nickel and/or titanium and/or nickel/titanium-based materials
- Recent pelvic venous thrombosis
- Recent myocardial infarction or bypass surgery within the last 30 days
- Shunt reversal with separated or significant right-to-left shunt

- Patients whose size or condition would cause the patient to be a poor candidate choice for cardiac catheterization (e.g. too small for echocardiographic imaging probe, catheter size, vasculature size, active infection, body weight < 8 kg)
- Patients with rim sizes of < 5 mm to the coronary sinus, inferior vena cava rim, an atrioventricular valve, or the right upper pulmonary vein

## 4. WARNINGS

The Figulla® Flex II ASD Occluder must be implanted exclusively by physicians trained in its use and are experienced with interventional transcatheter ASD closure techniques.

- Physicians who implant the Figulla® Flex II ASD Occluder must be able to recognize, assess and manage procedure-associated emergencies. On-site cardiac surgical support with corresponding personnel must be available.
- Patients allergic to nickel and/or titanium and/or nickel/titanium-based materials may suffer an allergic reaction to this device. Certain allergic reactions can be serious; patients shall be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty in breathing, or inflammation of the face or throat.
- The Figulla® Flex II ASD Occluder is not intended for multifenestrated ASDs.
- The use of improperly size-matched devices could seriously affect hemodynamics and optimal results. Before using this device, physicians shall carefully review the hemodynamic parameters as well as sizing information printed on the labels of the Figulla® Flex II ASD Occluder and its accessories (pusher and Occlutech® Delivery Set, ODS). Physicians shall also review the sizing and compatibility chart in Section 9 (see Table 3: Device Sizes and Recommended Delivery Systems) of this Instructions for Use (IFU) before starting a procedure. Careful consideration shall be given to ensure accurate size matching of accessories (pusher and ODS) with the corresponding device (i.e. review of device labels and color coding).
- Before using the Figulla® Flex II ASD Occluder, the physician shall carefully review Section 8 (Implantation Procedure) including relevant instructions therein on connecting the Figulla® Flex II ASD Occluder to the pusher. If a secure connection between the Figulla® Flex II ASD Occluder and pusher is not verified as described, disconnection of the device from its pusher may occur and the Figulla® Flex II ASD Occluder may embolize causing a life-threatening situation.
- The Figulla® Flex II ASD Occluder should not be used with delivery systems other than those recommended in Section 9 (see Table 3: Device Sizes and Recommended Delivery Systems) of this IFU.
- After deployment and release of the Figulla® Flex II ASD Occluder, complications such as device dislocation or embolization may occur as result of erroneous positioning or sizing of the device. These complications can present a life-threatening situation to the patient.
- Embolized Figulla® Flex II ASD Occluder must be retrieved. An emergency kit for the retrieval of the Figulla® Flex II ASD Occluder must be available in the catheterization laboratory during the procedure.
- The Figulla® Flex II ASD Occluder should only be released from its pusher after the physician has ascertained that the device is positioned correctly. This should be determined by performing fluoroscopy and/or Transesophageal (TEE) or Intracardiac Echocardiography (ICE) to visualize the Figulla® Flex II ASD Occluder and to confirm that the device is positioned properly.
- Do not release the device from the pusher if the Figulla® Flex II ASD Occluder does not conform to its original shape; if the device position is unstable, or if the device interferes with any adjacent cardiac structures such as the superior vena cava (SVC), pulmonary vein (PV), mitral valve (MV), coronary sinus (CS) or aorta (AO). In this case a correction of the position is required. In this case, please follow the guidance in Section 8.3: Notes for

Troubleshooting, “incorrectly positioned device”. If it is not possible to improve the position of the Figulla® Flex II ASD Occluder, the device must be removed and discarded.

- While still connected to the pusher, the Figulla® Flex II ASD Occluder can be retrieved or re-positioned using the recommended delivery system. When retrieving the device, do not move the Figulla® Flex II ASD Occluder through any intracardiac structure as this may cause damage to the intracardiac structure and can present a life-threatening situation to the patient.
- The Figulla® Flex II ASD Occluder must be used exclusively in accordance with these Instructions for Use and its implantation is to be carried out as described in this IFU.
- The physician shall inspect all packaging and labels of all devices before opening and follow the instructions. If the product box or sterile packaging is damaged in any manner, the Figulla® Flex II ASD Occluder shall be considered as unsterile and should not be used.
- The physician shall not use this device or any of its components if a seal appears to be or is broken (contents may not be sterile); if the label appears marked with text or symbols other than those on the label shown in this IFU or if the label is illegible, inappropriate, or absent.
- The physician shall not use this device or any of its components after the “use by” (expiration) date.

The Figulla® Flex II ASD Occluder is intended for single use only and this medical device is not suitable for re-sterilization. As soon as the sterile seal of the device is broken, it is contaminated. Re-use or re-sterilization may compromise the structural integrity of the devices, lead to device failure, and result in patient injury, illness or death.

- If, after inspection by the physician, the Figulla® Flex II ASD Occluder appears damaged or does not appear to function properly, the device is not suitable for implantation and should be discarded.
- Patients should be advised to avoid strenuous physical activity for a period of at least two weeks after device implantation.
- Patients with a rim < 5 mm from the aortic root may have a higher risk of erosion and if closed using a device, they will require closer follow up.
- Patients with rim sizes < 5 mm to the coronary sinus, inferior vena cava rim, an atrioventricular valve, or the right upper pulmonary vein may have a higher risk of device embolization and it is best to avoid doing such cases.

## 5. PREVENTATIVE MEASURES

### 5.1 Handling

Prior to introducing the Figulla® Flex II ASD Occluder in the patient, it is recommended to check whether the device reverts to its original shape once it has been pulled in and pushed out of the loader while immersed in a sterile physiological NaCl solution. If the Figulla® Flex II ASD Occluder does not revert to its intended shape, it is not suitable for implantation and shall be discarded (see details in Section 8: Implantation Procedure). The physician shall inspect the device prior to its use in patient and check the following.

- Pusher size is correct (match color codes);
- The Figulla® Flex II ASD Occluder can be properly connected to the pusher;
- The Figulla® Flex II ASD Occluder can be released from the pusher by operating the release mechanism;
- The Figulla® Flex II ASD Occluder is firmly attached to the pusher by gently pulling the device a few times;

- The Figulla® Flex II ASD Occluder reverts to its original shape once it has been pulled in and pushed out of the loader while immersed in a sterile physiological NaCl-solution

In case of any kind of inconsistency, the physician should select another device.

## 5.2 Device Size

A suitable Figulla® Flex II ASD Occluder is determined by ASD assessment (see Section 8.1: Measurement of Defect and Device Sizes) and the device-related sizing information in Section 9 (see Table 3: Device Sizes and Recommended Delivery Systems).

## 6. ADVERSE EVENTS

- Air embolism
- Allergic reactions
- Anesthesia reactions
- Apnea
- Arrhythmias
- AV-Fistula
- Bleeding (hemorrhage) requiring treatment
- Cardiac / vascular perforation
- Cardiac tamponade
- Death
- Embolization (peri & post procedural)
- Esophageal injury
- Femoral access complication
- Fever
- Hematoma
- Hypertension or Hypotension
- (Immediate) surgical interventions
- Infections, including endocarditis
- Pericardial effusions
- Pericarditis
- Post-pericardiotomy syndrome
- Pseudo aneurysm
- Pulmonary edema
- Seizure
- Stroke
- Tissue erosion
- Thrombus formation on device
- Thrombosis
- Valvular regurgitation

## 7. MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

Non-clinical testing and MRI simulations were performed to evaluate the entire family of the Figulla® Flex II ASD Occluder.

Non-clinical testing demonstrated that the entire family of this implant is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W / kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under defined scan conditions, the Figulla® Flex II ASD Occluder is expected to produce a maximum temperature rise of 3.1 °C after 15-minutes of continuous scanning (i.e., per pulse sequence).

### Artifact Information

In non-clinical testing, the image artifact caused by the Figulla® Flex II ASD Occluder extends approximately 5-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

## 8. IMPLANTATION PROCEDURE

The Figulla® Flex II ASD Occluder is intended for transcatheter delivery and techniques applied by physicians can vary. To minimize the risk of adverse events occurring, the following manufacturer's recommendations must be reviewed before implantation of this device.

### 8.1 Measurement of Defect and Device Sizing

**Warning:** The use of improperly size-matched devices could seriously affect hemodynamics and optimal results. Before using this device, physicians shall carefully review the hemodynamic parameters as well as sizing information printed on the labels of the Figulla® Flex II ASD Occluder and its accessories (pusher and ODS). Physicians shall also review the sizing and compatibility chart in Section 9 (see Table 3: Device Sizes and Recommended Delivery Systems) of this IFU before starting a procedure. Careful consideration shall be given to ensure accurate size matching of its accessories with the corresponding device (i.e. review of device labels and color-coding).

- Access the right femoral vein and perform a routine right heart catheterization
- Assess the secundum-type ASD simultaneously either by TEE or ICE resulting in a comprehensive echocardiographic study of all aspects of ASD anatomy, including location, size, presence of additional defects and adequacy of the Superior/Superior Vena Cava rim; Anterior-Superior/Aortic rim; Inferior/Inferior Vena Cava and Coronary Sinus rim; and Posterior rim)
- Select a guide wire (.035-.038”), cross the atrial septal defect and place it in the left upper pulmonary vein
- Select a multipurpose catheter (e.g. 5-French MP catheter), cross the atrial septal defect and place it in the left upper pulmonary vein
- Prepare the appropriate Sizing balloon according to manufacturer`s guideline
- Advance the sizing balloon over the guidewire into the left atrium and place it across the ASD under fluoroscopic and echocardiographic guidance
- Inflate the balloon until a clear defined waist is visible (TEE and fluoroscopic guidance) and measure it
- **Defect Sizing by using the “stop-flow-technique”:** The balloon is inflated with diluted contrast fluid until the left-to-right shunt ceases as observed by color-flow Doppler TEE/ICE.
- Perform precise ASD measurement by echocardiographic and/or fluoroscopic visualization of indentation made by ASD margins (best ECHO view: observing the balloon in its long axis)

**Warning:** To mitigate the risk of erosion, physicians must accurately evaluate the ASD. Implantation of the Figulla® Flex II ASD Occluder in patients with a deficient retro-aortic rim less than 5 mm (“naked aorta”), or a deficient superior rim significantly increases the risk of erosion.

### 8.2 Procedural Method of Implantation

*Required Accessories (not included in Figulla® Flex II ASD Occluder Set)*

- Occlutech® Delivery Set (ODS) containing:
  - Delivery sheath and dilator, loader and Y-hemostasis valve

- Guidewire
- Appropriate accessories to manage potential adverse events

### ***Recommended Imaging Tools***

- Fluoroscopy
- TEE or ICE (recommended before, during and after implantation. If TEE is used, the patient's esophageal anatomy must be adequate for the placement and manipulation of the TEE probe)

### ***Recommended Medication During Implantation Procedure***

- 48 hours' prior the procedure, patients are asked to take acetylsalicylic acid (3 – 5 mg/ kg per day).
- We recommended heparin to achieve Activated Clotting Time (ACT) > 200 s at the time of device deployment.
- Standard endocarditis prophylactic therapy should be administered during the procedure.

### ***Implantation***

The Figulla® Flex II ASD Occluder shall only be implanted by an experienced and trained physician, and in a specialized catheterization laboratory.


**Note:** The use of TEE or ICE is required before, during, and after implantation. If TEE is used, the patient's esophageal anatomy must be adequate for the placement and manipulation of the TEE probe.

1. Based on precise ASD measurement and device-related sizing information given in Section 9 (see Table 3: Devices Sizes and Recommended Delivery Systems), select an appropriate Figulla® Flex II ASD Occluder and the size-matched ODS noting further intracardiac structures as well.
2. Prepare the delivery sheath as described in the ODS IFU followed by its introduction together with the dilator into the right atrium. Sheath placement shall be performed under the guidance of imaging method, like TEE, fluoroscopy or ICE.
3. Once in place, the dilator shall be retracted by 10 mm before further moving the delivery sheath along the guidewire.
4. After positioning of delivery sheath in the left atrium, guidewire and dilator can be entirely removed.

**Note:** Do not abruptly remove the dilator from the left atrium as the risk for causing an air embolism can occur.

5. Remove the pusher and the Figulla® Flex II ASD Occluder from sterile packaging and check for device integrity and pusher functionality.
6. Place the Figulla® Flex II ASD Occluder in a bowl filled with a sterile physiological NaCl solution. Follow the pusher and ODS IFU and instructions given in Table 1 (step 6a to 6f), to connect the Figulla® Flex II ASD Occluder to the pusher.

**Table 1:** Device-Pusher-Connection

Steps	Flex – Pusher II (FP II)	Occlutech® Pistol Pusher (OPP)
<p><b>6a.</b> Insert the distal end (jaw-container) of the pusher through the Y-hemostasis valve into the loader (both provided with the ODS; please see ODS IFU for detailed device components and descriptions).</p>		
<p><b>6b.</b> Open the jaws of the applied pusher (please see pusher IFU for detailed device components and descriptions).</p>	<p>Turn the FP II lock mechanism (shown in red) to the ‘unlock’ position.</p>  <p>Push the handle at the proximal end of the FP II forward into distal direction (opens jaws).</p> 	<p>Please be sure that ‘unlock’ button is in the correct unlock position.</p>  <p>Pull the trigger located in the handle of OPP to control unlock position of handle and to release the jaws while pusher is in a straight position.</p> 
<p><b>6c.</b> Connect the open pusher jaws to the ball-connector located on the Figulla® Flex II ASD Occluder</p>		
<p><b>6d.</b> Retracting the jaws to connect the Figulla® Flex II ASD Occluder</p>	<p>Release the FP II handle.</p>	<p>Release the trigger mechanism of the OPP.</p>
<p><b>6e.</b> Securing the device-pusher connection while the pusher is in a straight position.</p> <p><b>Note:</b> A secure device-pusher connection is required to prevent any unintentional release of the device from the pusher.</p>	<p>Turn the FP II lock mechanism (shown in red) to its ‘locked’ position.</p> 	<p>Push the lock button of the OPP to its ‘locked’ position.</p> 

Steps	Flex – Pusher II (FP II)	Occlutech® Pistol Pusher (OPP)
6f. Figulla® Flex II ASD Occluder / Pusher assembly		

**Note:** The Figulla® Flex II ASD Occluder / pusher connection must be securely locked to prevent accidental release of the device from the pusher: If using the FP II, its lock mechanism (see Table 1: shown in red) must be screwed tightly. If using the OPP, check that the lock mechanism is in its ‘locked’ position.

7. The Figulla® Flex II ASD Occluder is now attached to the pusher and shall be pulled into the loader while both, device and loader (distal end) are immersed in a sterile physiological NaCl-solution.
8. The following steps constitute a potential risk of introducing air into the delivery system which may lead to an air embolism. So it is important to ensure that all air-bubbles are removed from the loader and the delivery sheath by physician as described:
  - After the Figulla® Flex II ASD Occluder has been completely retracted into the loader, the loader / device assembly must be thoroughly flushed via the side arm of Y-valve using the sterile physiological NaCl-solution in order to remove any air-bubbles.
  - Before connecting the loader / device assembly to the delivery sheath, it is necessary to ensure that some blood leaks from the delivery sheath to remove any air bubbles (back bleeding).
9. Once the loader / device assembly is firmly connected to the delivery sheath, the Figulla® Flex II ASD Occluder is advanced from the loader into the delivery sheath by pushing the pusher up to the distal end of the delivery sheath (positioned in the left atrium).
10. Once in place, the left atrial (LA) disc and the waist of the Figulla® Flex II ASD Occluder are opened by carefully advancing the pusher towards the distal end of the delivery sheath. After the LA disc has been successfully deployed in the left atrium, the device is held in its position while slightly pulling it against the atrial septum by carefully retracting the delivery sheath. An elastic resistance with pulse synchronized movements and TEE or ICE imaging should be used to confirm proper positioning of LA disc at this time.
11. Once confirmed, the proximal disc of the Figulla® Flex II ASD Occluder is opened in the right atrium, by gently retracting the delivery sheath until the device has been completely deployed.
12. Correct positioning of the Figulla® Flex II ASD Occluder must be ensured by injecting contrast agent through the delivery sheath which should also be used to appropriately visualize any residual shunt or a sub-optimally fitting product.
13. If the Figulla® Flex II ASD Occluder does not adequately fit upon deployment, a correction of its position must be considered (see Section 8.3: Notes for Troubleshooting). If repositioning of the device cannot be achieved, or a significant shunt remains, the Figulla® Flex II ASD Occluder should be withdrawn back into the delivery sheath, removed and discarded.
14. If the optimal fit of the Figulla® Flex II ASD Occluder has been established and confirmed by using TEE or ICE, a gentle pull-and-push maneuver (tug test) can be considered to verify the stable positioning of the device.

15. If the Figulla® Flex II ASD Occluder passes the tug test, the device can be released from the pusher, otherwise it needs to be retrieved and repositioned (see Section 8.3: Notes for Troubleshooting).
16. During the Figulla® Flex II ASD Occluder release from the pusher (see Table 2, step 16a and b), the physician should hold the device in place with the aid of the ODS sheath. The distal tip of this sheath should be kept close to the proximal end of the device to ensure an easy retrieval in case of device dislocation.

**Table 2:** Releasing the device from the pusher

Steps	Flex – Pusher II (FP II)	Occlutech® Pistol Pusher (OPP)
16a. Unlocking the Device-pusher-connection	Turn the lock mechanism (shown in red) of FP II to the ‘unlock’ position. 	Move the unlock button to the handle of the OPP to its ‘unlocked’ position. 
16b. Open the jaws and disconnect / release the Figulla® Flex II ASD Occluder.	Push the handle at the proximal FP II end forward. 	Slowly pull the trigger mechanism located in the handle at the proximal OPP end. 

**Note:** The physician should only disconnect the Figulla® Flex II ASD Occluder from the pusher once proper positioning of the device has been confirmed using echocardiographic imaging.

17. Now, slowly retract the pusher into the ODS sheath, entirely remove both and discarded them (see Section 11: Disposal).

### 8.3 Notes for Troubleshooting

#### *Incorrectly positioned device*

If the Figulla® Flex II ASD Occluder is not properly positioned after opening both discs, a correction of the position is required. To do so, the device must be completely withdrawn into the delivery sheath. Thereafter, another attempt to place and position the Figulla® Flex II ASD Occluder can be made. If positioning of the device cannot be improved, the Figulla® Flex II ASD Occluder should be withdrawn into the delivery sheath, removed and discarded.

### *Misconfigured device*

If the Figulla® Flex II ASD Occluder does not develop its intended shape during positioning, the device must be withdrawn completely into the delivery sheath by retracting the pusher. The Figulla® Flex II ASD Occluder should be removed and discarded.

## **8.4 Report of Serious Adverse Events**

All serious adverse events which are observed in relation to any Occlutech® devices have to be reported to the manufacturer and the corresponding national authority by physician or a delegated person (Regulation 2017/745 on Medical Devices, Chapter III of Annex I).

## **8.5 Recommended Follow-up Examinations**

### 1st day after implantation

- X-ray screening for status check
- TEE, TTE or ICE (recommended; preclusion of left atrial thrombus or thrombus adhering to the device, preclusion of residual shunt)
- Electrocardiography (ECG)
- Laboratory values

### 30th day after implantation

- TEE, TTE or ICE ((preclusion of left atrial thrombus or thrombus adhering to the device, preclusion of residual shunt)
- Clinical-neurological examination
- ECG for rhythm documentation

### 60th day after implantation

- TEE, TTE or ICE (preclusion of left atrial thrombus or thrombus adhering to the device, preclusion of residual shunt)

### 180th day after implantation

- TEE, TTE or ICE (recommended; preclusion of left atrial thrombus or thrombus adhering to the device, preclusion of residual shunt)
- ECG.

### 360th day after implantation

- TEE, TTE or ICE (preclusion of left atrial thrombus or thrombus adhering to the device, preclusion of residual shunt)
- ECG.


## **8.6 Recommended Medication after Implantation**

- Standard endocarditis prophylactic therapy should be administered for a minimum of 6 months' post implantation.
- Antithrombotic therapy is indicated regarding patient conditions. Clopidogrel (75 mg daily) for 2 - 3 months and acetylsalicylic acid (81 – 325 mg orally once daily) for 6 months' post-implantation is recommended.

## 9. DEVICE SIZES AND RECOMMENDED DELIVERY SYSTEMS

### Compatible systems:

- Pusher: Flex – Pusher II (FP II) or Occlutech® Pistol Pusher (OPP)
- Delivery Set: Occlutech® Delivery Set with delivery sheath, dilator, loader and Y-hemostasis valve

Ø LA Disc: Left Atrial Disc Diameter,  
 Ø RA Disc: Right Atrial Disc Diameter  
 Ø Waist: Max. Waist Diameter

Figure 2: Representative picture of Figulla® Flex II ASD Occluder

Table 3: Device Sizes and Recommended Delivery Systems

### Figulla® Flex II ASD Occluder

Figulla® Flex II ASD Article no.	Occlutech® Delivery Set Size [F]	Flex – Pusher II Article no.	Occlutech® Pistol Pusher Article no.	Ø Waist [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	AS Defect Size (D) [mm]
29ASD04	7 F	51FP060 (green)	55PP090 (green)	4	11	9	$D \leq 4$
29ASD05	7 F	51FP060 (green)	55PP090 (green)	5	14	11	$4 < D \leq 5$
29ASD06	7 F	51FP100 (light blue)	55PP125 (light blue)	6	16.5	12.5	$5 < D \leq 6$
29ASD07	7 F	51FP100 (light blue)	55PP125 (light blue)	7.5	18	14	$6 < D \leq 7.5$
29ASD09	7 F	51FP100 (light blue)	55PP125 (light blue)	9	20.5	16.5	$7.5 < D \leq 9$
29ASD10	7 F	51FP100 (light blue)	55PP125 (light blue)	10.5	22	18	$9 < D \leq 10.5$
29ASD12	9 F	51FP120 (yellow)	55PP165 (yellow)	12	27	23	$10.5 < D \leq 12$
29ASD13	9 F	51FP120 (yellow)	55PP165 (yellow)	13.5	28.5	24.5	$12 < D \leq 13.5$
29ASD15	9 F	51FP120 (yellow)	55PP165 (yellow)	15	30	26	$12 < D \leq 15$
29ASD16	9 F	51FP120 (yellow)	55PP165 (yellow)	16.5	31.5	27.5	$15 < D \leq 16.5$
29ASD18	9 F	51FP120 (yellow)	55PP165 (yellow)	18	33	29	$15 < D \leq 18$
29ASD19	10 F	51FP120 (yellow)	55PP165 (yellow)	19.5	34.5	30.5	$16.5 < D \leq 19.5$

## *Figulla® Flex II ASD Occluder*

<b>Figulla® Flex II ASD Article no.</b>	<b>Occlutech® Delivery Set Size [F]</b>	<b>Flex – Pusher II Article no.</b>	<b>Occlutech® Pistol Pusher Article no.</b>	<b>Ø Waist [mm]</b>	<b>Ø LA Disc [mm]</b>	<b>Ø RA Disc [mm]</b>	<b>AS Defect Size (D) [mm]</b>
29ASD21	11 F	51FP150 (purple)	55PP185 (purple)	21	36	32	18 < D ≤ 21
29ASD24	11 F	51FP150 (purple)	55PP185 (purple)	24	39	35	21 < D ≤ 24
29ASD27	12 F	51FP160 (blue)	55PP210 (blue)	27	42	38	24 < D ≤ 27
29ASD30	12 F	51FP160 (blue)	55PP210 (blue)	30	45	41	27 < D ≤ 30
29ASD33	12 F	51FP160 (blue)	55PP210 (blue)	33	48	43	30 < D ≤ 33
29ASD36	12 F	51FP160 (blue)	55PP210 (blue)	36	52	46	33 < D ≤ 36
29ASD39	12 F	51FP160 (blue)	55PP210 (blue)	39	54	49	36 < D ≤ 39
29ASD40	12 F	51FP160 (blue)	55PP210 (blue)	40	55	50	39 < D ≤ 40

## 10. STORAGE CONDITIONS


Temperature: 5 to 30 °C

## 11. DISPOSAL

All Occlutech® devices have to be disposed of properly according to the applicable regulations, guidelines or instructions on waste disposal given by the individual clinic.

**Notes:** The empty outer and sterile packaging and IFUs can be separately put in the corresponding collection containers for recycling purposes.

## 12. LABELLING

	CAUTION		MANUFACTURER ADDRESS
	DO NOT USE IF PACKAGE IS DAMAGED.		MANUFACTURING DATE
	DO NOT RESTERILIZE.		CONSULT INSTRUCTION FOR USE
	SINGLE USE ONLY.		SERIAL NUMBER (UNIQUE PRODUCT ID-NUMBER )
	STERILIZED WITH ETHYLENE OXIDE.		ITEM NUMBER / ORDER NUMBER
	CE MARK. PRODUCT MEETS THE REQUIREMENTS ACCORDING TO 93/42/EEC		STORE DRY & AWAY FROM SUNLIGHT
	USE BY DATE. / EXPIRE DATE.		TEMPERATURE LIMIT FOR (LONG TERM) STORAGE
	DIAMETER [mm]	<b>Quantity</b> <b>1</b>	QUANTITY

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**Occlutech GmbH